## Cell Cap Removal Tool



**Brand:** 

**Product Code: 170-33** 

Phone: 832-320-7300 - Email: sales@ofite.com

**Availability:** 16 **Lead Time:** 60

## **Description**

The Cell Cap Removal Tool (#170-33) provides a convenient method for removing a cell cap from a HTHP filter press cell. After high temperature testing, cell caps are often difficult to remove due to expansion and o-ring friction. This tool screws into the valve stem port and pivots against the side wall of the cell. Simply pushing down on the lever pulls the cap directly up and out of the cell body.